EN Dymedix FastTrack Interface Kit User Guide

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
EN

CAUTION STATEMENTS

In the event of skin reaction, discontinue use.

This product is for diagnostic purposes only and is not to be used in life supporting or life sustaining situations.

Any serious incident that has occurred in relation to this device should be reported to Dymedix Diagnostics and the competent authority of the Member States in which the user and/or patient is established.

Symbol	Definition of Symbol
	Latex Free
	Do not re-use
	Federal Law (USA) restricts this device to sale by or on the order of a physician.
	Medical Device
	Consult instructions for use.
	Date of Manufacture Expressed as: YYYYMMDD
	Manufactured by

FEATURES AND INTENDED USE

The Dymedix FastTrack Interface Cable connects the FastTrack Effort Pack sensors to the PSG headbox during sleep disorder studies. Features and specifications of this device are subject to change without notice. This Interface Kit must be used with a Dymedix FastTrack Effort Pack.

INSTRUMENT CONNECTIONS & SETTINGS

Refer to the polysomnogram operator’s manual for direction in determining proper instrument connection and settings. Response is affected by variables such as sensitivity settings, high and low frequency filter settings, sensor placement, patient movement, battery life etc. Incorrect connection could result in erratic or no signal.Adjusting the sensitivity (gain) up or down is usually required.

PRECAUTIONS

- Store FastTrack Interface Kit cables at room temperature.
- Do not sterilize Interface Cables.
- Use all available recorded parameters in making diagnosis.

PATIENT APPLICATION

- Connect the FastTrack effort sensors to the two receiving connections on the Interface Cable matching the red/black color coded ends.
- Connect the Interface Cable to the headbox of the PSG or HST device. Match colors of keyhole connectors to the appropriate effort inputs on the headbox. If using interface cables with non-keyhole connectors, connect red/black DIN connectors for correct polarity.
- The Interface Cable must be connected to an AC-coupled channel on the recording instrument. Proceed with polysomnography as outlined by clinician procedure.
- If there is no signal output from the sensor and it is properly placed on the patient discontinue use.

REMOVAL INSTRUCTIONS

- Disconnect Interface Cable leads from recording equipment.
- The Interface Cable may be cleaned with a non-corrosive hospital-approved disinfectant.
- Do not immerse in liquids.**

WARRANTY


Dymedix warrants its sensors to be free from defects in materials and workmanship and to perform satisfactorily for One (1) Year from the date of purchase. Should the Interface Cable fail during the warranty period, the sole liability of Dymedix is limited to replacement of the product with no charge for parts or labor. Under no circumstances shall Dymedix be liable for any consequential damages arising from the use of, or the inability to use this product. This warranty is rendered void for conditions resulting from: damage, marginal performance or malfunctions caused by: misuse, abuse, neglect, improper amplification settings, improper line voltage, power fluctuations, or any adverse environmental conditions, tampering, unauthorized modifications, adjustments or repairs to the product or its accessories. This warranty is in lieu of all other warranties, expressed or implied and is extended only to the original purchaser.

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