EN

Perfect Fit™ Adult Effort Belt Kit

FR

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
EN

CAUTION STATEMENTS

In the event of skin reaction, discontinue use.

This product is for diagnostic purposes only and is not to be used in life supporting or life sustaining situations.

Any serious incident that has occurred in relation to this device should be reported to Dymedix Diagnostics and the competent authority of the Member States in which the user and/or patient is established.

Symbol	Definition of Symbol
	Latex Free
	Do not re-use
	Federal Law (USA) restricts this device to sale by or on the order of a physician.
	Medical Device
	Consult instructions for use.
	Date of Manufacture Expressed as: YYYYMMDD
	Manufactured by

FEATURES AND INTENDED USE

The Dymedix Perfect Fit™ respiratory belt sensors have been designed to detect respiratory effort signals during sleep studies for recording by a physiological recorder. These sensors produce electrical signals when respiratory effort occurs. The signal size can be adjusted by changing the settings on the recording equipment. The use of a Dymedix EFM (Electronic Filtration Module) is required with this product. Features and specifications of this device are subject to change without notice. The Dymedix respiratory effort belt sensor is used with existing sleep study-recording devices in support of diagnostic recording of respiratory effort. Discontinue use of device when signal amplitude falls below recordable levels or the unit is damaged.

INSTRUMENT CONNECTIONS & SETTINGS

Refer to the polysomnograph operator's manual for direction in determining proper instrument connection and settings. Response is affected by variables such as sensor placement, sleep position, ambient room temperature, sensitivity settings, frequency filter settings, EFM battery life and system characteristics.

PRECAUTIONS

1. Store sensor at room temperature.
2. Do not sterilize sensor.
3. Connect the sensor to electrically isolated equipment.
4. Incorrect connections could result in erratic or no signal.
5. Use all available recorded parameters in making diagnosis.
6. Do not apply electrode to broken or damaged skin.

PATIENT APPLICATION

1. Dymedix recommends using Perfect Fit™ Disposable Effort Belts on any patients with suspected Infectious Disease.
2. Wrap the proper belts around the chest and abdominal areas with the Perfect Fit™ logo facing away from patient. The belts should be worn over night clothes.
3. Attach an appropriately sized Perfect Fit™ extender strap to each sensor and adjust to proper tension. **DO NOT OVER-TIGHTEN BELT.** Belts should be snug but still allow a hand to slide in under the belt.
4. Insert yellow/black abdomen belt connectors into yellow/black EFM cable end. Repeat for blue/black chest connectors.
5. Connect the EFM to the headbox per EFM instructions.
6. For proper strain relief, the use of stress loops is recommended.
7. Proceed with polysomnography test as outlined by clinician procedure.
8. **Note: BE SURE PATIENT IS PROPERLY GROUNDED.**
9. If there is no signal output from the sensor and it is properly placed on the patient discontinue use.
- Note:** Observe that an upward waveform deflection occurs on inhalation and a downward waveform deflection occurs on exhalation. Reverse the input connections, if necessary, to obtain proper waveform deflection.

SENSOR REMOVAL INSTRUCTIONS

1. Disconnect EFM cable from the belt.
2. Remove belt from patient.

Caution: Firmly hold connectors to remove. Pulling wires will cause damage that is not covered under warranty.

CLEANING INSTRUCTIONS

1. Clean the 8" sensor and wires with non-corrosive hospital approved fabric cleaning agent and let dry. Do not immerse the sensor. Sensor must be completely dry before using on patient.
2. Clean the reusable belt straps with a non-corrosive hospital approved fabric laundry detergent. Completely rinse the reusable belt strap in water.
3. Machine-washing straps on a gentle cycle is also permitted.
4. Allow the reusable belt straps to dry. Belt straps must be completely dry before using on a patient.
5. Do not use bleach to clean.

WARRANTY

Dymedix warrants its sensors to be free from defects in materials and workmanship and to perform satisfactorily for One (1) Year from the date of purchase. Should the sensor fail, the sole liability of Dymedix is limited to replacement of the product with no charge for parts or labor. Under no circumstances shall Dymedix be liable for any consequential damages arising from the use of, or the inability to use this product. This warranty is rendered void for conditions resulting from: damage, marginal performance or malfunctions caused by: reuse, misuse, abuse, neglect, improper amplification settings, improper line voltage, power fluctuations, or any adverse environmental conditions, tampering, unauthorized modifications, adjustments or repairs to the product or its accessories. This warranty is in lieu of all other warranties, expressed or implied and is extended only to the original purchaser.

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