DyMedix[®]

DIAGNOSTICS

EN Snore Sensor, Disposable

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DA Snorkesensor til engangsbrug

NL Snurksensor, wegwerp

DE Einweg-Schnarchsenssor

IT Sensore del russamento, monouso

NO Snorkesensor for engangsbruk

PT Sensor de Ronco, Descartável

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EN Snore Sensor, Disposable

86-41021-0010 Rev F

EN Connect the sensor wires to the EFM cable. Connect the EFM cable to the input connections on the recording instrument head box.

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EN

CAUTION STATEMENTS

In the event of skin reaction, discontinue use. Incorrect placement or connections could result in erratic or no signal.

This product is for diagnostic purposes only and is not to be used in life supporting or life sustaining situations.

Any serious incident that has occurred in relation to this device should be reported to Dyemedix Diagnostics and the competent authority of the Member States in which the user and/or patient is established.

Symbol	Definition of Symbol
	Latex Free
	Do not re-use
	Federal Law (USA) restricts this device to sale by or on the order of a physician.
	Medical Device
	Consult instructions for use.
	Date of Manufacture Expressed as: YYYYMMDD
	Manufactured by

FEATURES AND INTENDED USE

The Dyemedix Snore sensor detects snoring sound signals during sleep disorder studies. This sensor is intended to produce a small electrical signal to existing recording equipment. The signal size can be adjusted by changing the settings on the recording equipment. This sensor must be used with a Dyemedix Electronic Filtration Module (EFM)

INSTRUMENT CONNECTIONS & SETTINGS

Refer to the polysomnograph manual for proper instrument connection and settings. Response is affected by sensor placement, sleep position, ambient room temperature, sensitivity settings, frequency filter settings, EFM battery life and system characteristics.

PRECAUTIONS FOR SENSOR

- Store sensor at room temperature.
- Do not sterilize or reuse sensor.
- Connect the sensor to electrically isolated equipment.
- Always disconnect the sensor from the EFM cable when data is not being collected.
- Use all available recorded parameters in making diagnosis.
- If there is no signal output from the sensor and it is properly placed on the patient discontinue use.
- If the unit is damaged or the adhesive will not stick discontinue use and discard.

PATIENT INSTALLATION

- Remove the release liner from the back of the sensor to expose the adhesive surface.
- Sensor should fit as shown below, as the patient is making snoring sounds, gently place your fingers on the fron of the patient's neck. Apply the sensor to the clean and unbroken skin over the patient's larynx where you feel the most vibration.
- Connect the sensor wires to the EFM cable. Connect the EFM cable to the input connections on the recording instrument head box.
- The sensor must be connected to AC-coupled inputs on the recording instrument.
- BE SURE PATIENT IS PROPERLY GROUNDED.**
- Proceed with polysomnography as outlined by clinician procedure.
- Do not apply electrode to broken or damaged skin.

SENSOR REMOVAL INSTRUCTIONS

- Disconnect sensor wires from EFM cable.
- Remove sensor from patient.
- Dispose of sensor as normal waste. Do not reuse.

WARRANTY

Dyemedix warrants its sensors to be free from defects in materials and workmanship and to perform satisfactorily for One (1) Year from the date of purchase. Should the sensor fail, the sole liability of Dyemedix is limited to replacement of the product with no charge for parts or labor. Under no circumstances shall Dyemedix be liable for any consequential damages arising from the use of, or the inability to use this product. This warranty is rendered void for conditions resulting from: damage, marginal performance or malfunctions caused by: reuse, misuse, abuse, neglect, improper amplification settings, improper line voltage, power fluctuations, or any adverse environmental conditions, tampering, unauthorized modifications, adjustments or repairs to the product or its accessories. This warranty is in lieu of all other warranties, expressed or implied and is extended only to the original purchaser.

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MedEnvoy Global B.V. Prinses Margrietplantsoen 33 Suite 123 2595AM The Hague The Netherlands