SnuggleFit™ Infant Effort Belt System Reusable Effort Module

EN Dymedix SnuggleFit Infant Effort Belt

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86-41600-0110 Rev B © 2019 Copyright. All rights reserved.

SnuggleFit™ Infant Effort Belt Dymedix is a registered trademark of Dymedix Diagnostics Inc.

EN

CAUTION STATEMENTS

This product is for diagnostic purposes only and is not to be used in life supporting or life sustaining situations.

Any serious incident that has occurred in relation to this device should be reported to Dymedix Diagnostics and the competent authority of the Member States in which the user and/or patient is established.

Symbol	Definition of Symbol
	Latex Free
	Do not re-use
	Federal Law (USA) restricts this device to sale by or on the order of a physician.
	Medical Device
	Consult instructions for use.
	Date of Manufacture Expressed as: YYYYMMDD
	Manufactured by

FEATURES AND INTENDED USE

The Dymedix SnuggleFit™ Infant Effort sensors have been designed to detect respiratory effort signals during sleep studies for recording by a physiological recorder. These sensors produce electrical signals when respiratory effort occurs. The signal size can be adjusted by changing the settings on the recording equipment. The SnuggleFit™ Infant Effort Belt Module must be used with the Dymedix disposable button-hole effort straps. The use of a Dymedix SnuggleFit™ Interface Kit is required with this product. Features and specifications of this device are subject to change without notice.

INSTRUMENT CONNECTIONS & SETTINGS

Refer to the PSG operator’s manual for direction in determining proper connections and settings. Response maybe affected by variables such as sensor placement, sleep position, ambient room temperature, sensitivity setting, and frequency filter settings.

PRECAUTIONS

- Store sensor at room temperature.
- Do not sterilize sensor.
- Connect the sensor to electrically isolated equipment.
- Use all available recorded parameters in making diagnosis.
- Do not apply electrode to broken or damaged skin.

PATIENT APPLICATION

- Place an optional animal sticker on each SnuggleFit™ Effort Sensor.
- Wrap the individual effort straps around the chest and abdominal areas. The belts should be worn over sleeping clothes.
- Determine appropriate length and connect the ends of the strap using the button connector. Do not over tighten. Straps should be snug but still allow a hand to slide in under the belt.
- Any remaining excess strap material can be cut with a scissor and discarded. However, it is recommended that a 6 inch length of material be left for any future adjustments that might be necessary during the test.
- Place an optional SnuggleFit™ Comfort Wrap around each sensor for those patients who may be pressure sensitive.
- Using the SnuggleFit™ Interface Kit, connect both straps to the interface cable inputs, and plug the interface cable output jacks into the recording device. (Be sure to match color coded inputs and outputs to ensure correct waveform polarity).
- If there is no signal output from the sensor and it is properly placed on the patient discontinue use.

SENSOR REMOVAL INSTRUCTIONS

- Remove both straps.
- Disconnect the Interface Cables from the straps.
- Remove the reusable Interface Cable from the PSG. **Caution: Firmly hold connectors when inserting or removing cables. Pulling wires will cause damage that is not covered by warranty.**
- Dispose of sensor straps as normal waste. Do not reuse.

WARRANTY

Dymedix warrants its sensors to be free from defects in materials and workmanship and to perform satisfactorily for One (1) Year from the date of purchase. Should the Reusable Effort Module fail during the warranty period, the sole liability of Dymedix is limited to replacement of the product with no charge for parts or labor. Under no circumstances shall Dymedix be liable for any consequential damages arising from the use of, or the inability to use this product. This warranty is rendered void for conditions resulting from: damage, marginal performance or malfunctions caused by: misuse, abuse, neglect, improper amplification settings, improper line voltage, power fluctuations, or any adverse environmental conditions, tampering, unauthorized modifications, adjustments or repairs to the product or its accessories. This warranty is in lieu of all other warranties, expressed or implied and is extended only to the original purchaser.

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