

ALWAYS INNOVATING™
DyMedix®
DIAGNOSTICS
DYMEDIX.COM • 888.212.1100

SnuggleFit™
INFANT RESPIRATORY BELT

TriplePlay™
AIRFLOW SENSOR

PerfectFit™
EFFORT BELT 3D

iFM
INFANT FILTRATION MODULE

FM1
AIRFLOW FILTRATION MODULE

FM2
AIRFLOW FILTRATION MODULE

PerfectFit™
EFFORT BELT 2

FastTrack
EFFORT PACK

COMPATIBLE WITH ALL PSG AND EEG SYSTEMS | AASM AND JACHO COMPLIANT

MK-0089 APR 2020

DYMEDIX.COM | (888) 212-1100



Lullaby™
PEDIATRIC SENSORS
by **DyMedix**®
DIAGNOSTICS

The most advanced and innovative family of diagnostic sleep sensors made especially for pediatric patients.



PEDS ARE NOT LITTLE ADULTS.

Dymedix Diagnostics designs sleep sensors with comfort and accuracy in mind for infants and children, our most sensitive patients. Using PVDF technology, respiratory effort, movement and airflow monitoring has never been more accurate and easier for technologists and patients. Our Lullaby™ family of pediatric sensors are available in both disposable and reusable options.

AIRFLOW

TriplePlay™
AIRFLOW SENSOR

As part of the Lullaby™ sensor family, Dymedix Diagnostics TriplePlay® airflow sensors respond to both temperature and pressure changes providing an extremely sensitive airflow monitoring system for neonates, infants and children.

DISPOSABLE AIRFLOW SENSORS

Available in Neonate, Infant and Pediatric



TRIPLECHIN™ EMG PATCH

REUSABLE AIRFLOW SENSORS

Available in Pediatric and Infant



Our FM interface cables incorporate both passive and active filtration modules to help isolate specific airflow waveforms of interest.



RESPIRATORY EFFORT

Plug and Play, Throw it Away (or use it again)! The Lullaby™ family of respiratory effort systems provides pediatric sleep clinicians with reusable, semi-disposable or completely disposable options, all designed with the ease of use, slim profile and enhanced sensitivity required for quality pediatric respiratory effort monitoring.

DISPOSABLE RESPIRATORY EFFORT

SnuggleFit™
INFANT RESPIRATORY BELT



PerfectFit™
EFFORT BELT 3D



FastTrack™
EFFORT PACK



REUSABLE RESPIRATORY EFFORT

PerfectFit™
EFFORT BELT 2

