

We've Got You.



DISTRIBUTORS

Nihon Kohden America
1-800-325-0283
www.us.nihonkohden.com

MVAP Medical Supplies Inc.
(877) 735-6827
mvapmed.com

Electra-Medical Corp.
(800) 678-4856
electramed.com

Carrier Medical
(800) 960-1434
carriermedical.com

Compumedics USA Inc.
(877) 717-3975
compumedics.com.au

DRS4DRS
(866) 374-4377
drs4drs.com

Cadwell Industries Inc.
(800) 245-3001
cadwell.com

KEGO USA Inc.
(800) 600-1390
kegousa.com

Natus Medical Inc.
(800) 960-1434
natus.com

Bird Healthcare Pty Ltd.
(300) 365-561
birdhealthcare.com

Dy^Mmedix[®]

DIAGNOSTICS

DISPOSABLES

DISPOSABLE SLEEP SENSORS TO KEEP PATIENTS AND STAFF SAFE

ALWAYS INNOVATING™
Dy^Mmedix[®]
DIAGNOSTICS
DYMEDIX.COM • 888.212.1100

MK-0129

COMPATIBLE WITH ALL PSG & EEG SYSTEMS | AASM & JCAHO COMPLIANT

SIMPLE. SENSITIVE. DISPOSABLE.

Dymedix Diagnostics carries a full line of disposable respiratory effort, airflow, EEG, EMG and movement sensor products for the sleep and neurodiagnostic community. Dymedix Diagnostic's exclusive line of disposable sleep sensors are based on our patented technology called polyvinylidene fluoride (PVDF) which delivers very fast, and extremely reliable signals.


DISPOSABLE EMG & MOVEMENT SENSORS



ACCUSNORE™ is the ONLY disposable snore sensor available. Using PVDF technology the AccuSnore™ sensor provides superior sensitivity compared to conventional piezo snore sensors. Adhesive backing allows for multiple placements during testing to ensure accurate capture of snoring activity.



MOTIONWAVE™ DUAL EMG PATCH is the only completely disposable pre-measured dual EMG electrode patch. No snaps and no lead wires to clean. Use it once and discard the entire sensor, wires and all.



TRIPLECHIN™ EMG PATCH A completely disposable chin EMG electrode patch. Three electrodes perfectly spaced with an adhesive backing for easy application and reliable EMG signals. Available in both adult and pediatric sizes.

DISPOSABLE RESPIRATORY EFFORT BELTS



FAST TRACK™ EFFORT PACK is the first and only completely disposable respiratory effort solution. Each FastTrack™ Effort Pack includes single use Chest and Abdomen effort belts incorporating our cut-to-size Button-Hole™ strap system. Plug and Play, Throw it Away.



PERFECTFIT 3D™ EFFORT SYSTEM is an affordable respiratory effort system that incorporates PVDF effort modules with Button-Hole™ connectors to where disposable effort belt material is connected. PerfectFit 3D™ is a perfect semi-disposable solution for PSG and HST respiratory effort applications. Half the cost of RIP.


SNUGGLEFIT™ INFANT EFFORT SYSTEM is a disposable respiratory effort system designed specifically for infants. The SnuggleFit™ responds to both strain and pressure making it extremely sensitive, and incorporates our disposable cut-to-size Button-Hole™ disposable belt material.

DISPOSABLE AIRFLOW SENSORS



TRIPLEPLAY™ PVDF AIRFLOW SENSORS Extreme sensitivity and response time. Responds to changes in both temperature and pressure providing accurate detection of Apnea, Hypopnea, UARS and Snore, all with one sensor. Designed for use with a TriplePlay FM airflow filtration module.

DISPOSABLE EEG ELECTRODES



PRECISION GOLD™ DISPOSABLE EEG ELECTRODES Disposable gold-plated EEG cup electrodes at a very affordable price. Designed with a flat rim and rounded edges for exceptional patient comfort, especially during extended monitoring. Includes a flexible and durable hub connector with strain relief. Available in 4 different lengths.



CompletePlus

DISPOSABLE PSG SETUP KIT

Get the complete collection of Dymedix disposable sensors in one grab-and-go PSG setup kit. CompletePlus kits are available individually or in 20 or 40-kit boxes.

SEE THE ENTIRE DISPOSABLES COLLECTION AT DYMEDIX.COM.